

ZBTB25 antibody



Catalog Number: 116912

Product name

ZBTB25 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZBTB25 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse spleen tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of ZBTB25 recombinant protein (Accession Number: NM_001354686). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116912(ZBTB25 Antibody) at dilution of 1:600



Immunofluorescent analysis of HepG2 cells using Catalog No:116912(ZBTB25 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)