

Anti-PAX6 Rabbit antibody



Catalog Number: 168912

Product name

Anti-PAX6 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

PAX6 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human PAX6

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

1:1000

1:20-1:500

1:50

Validations



Western blot detection of PAX6 in K562, HeLa cell lysates using PAX6 Rabbit mAb(1:1000 diluted). Predicted band size: 47kDa. Observed band size: 47kDa.