

Anti-CEACAM1 antibody



Catalog Number: 141212

Product name

Anti-CEACAM1 antibody

Immunogen

[Rat CEACAM1/CD66a \(His Tag\) recombinant protein](#)

Specificity

Rat CEACAM-1/CD66a

Antibody description

Mouse Monoclonal to Rat CEACAM1

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Rat CEACAM-1/CD66a and conjugated with PE under optimum conditions, the unreacted PE was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from

prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

FCM

Dilutions

FCM

Validations



Rat CEACAM-1/CD66a Flow Cytometry (FC) 15888