

Anti-IRF5 Rabbit antibody



Catalog Number: 167212

Product name

Anti-IRF5 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

IRF5 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human IRF5

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP, FC

Dilutions

WB: 1/500-1/1000


IHC: 1/50-1/200

ICC/IF: 1/100-1/500

IP: 1/100

FC: 1/50

Validations



Western blot analysis of IRF5 expression in THP-1 cell lysate.