

# Anti-Ikaros (1A12) Mouse antibody



Catalog Number: 168212

## Product name

WB, IP

Anti-Ikaros (1A12) Mouse antibody

## Dilutions

## Specificity

WB: 1/1000

Human

## Antibody description

Ikaros (1A12) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human Ikaros(C-terminus) protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality


Monoclonal

## Ig Type


Mouse IgG1

## Applications

## Validations



Immunoprecipitation analysis of HeLa cell lysate using Ikaros (C-terminus) mouse mAb.



Western blot detection of Ikaros(C-terminus) in Ramos and Jurkat cell lysates using Ikaros(C-terminus) mouse mAb (1:1000 diluted). Predicted band size: 58KDa. Observed band size: 58KDa.