

ABCG8 Polyclonal Antibody



Catalog Number: 163612

Product name

ABCG8 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ABCG8

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human ABCG8 (NP_071882.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations



Western blot - ABCG8 Polyclonal Antibody

Western blot analysis of extracts of A-549 cells, using ABCG8 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - ABCG8 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using ABCG8 antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.