

SSTR1 antibody

Catalog Number: 115612

Product name

SSTR1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

SSTR1 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, BxPC-3 cells, mouse brain tissue, SH-SY5Y cells, mouse small intestine tissue. Observed molecular weight by Western-blot: 53kd,63kd

Preparation

This antibody was obtained by immunization of SSTR1 recombinant protein (Accession Number: NM_001049). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:115612(SSTR1 antibody) at dilution of 1:600