

Anti-LOC4331832 antibody



Catalog Number: 176612

Product name

Anti-LOC4331832 antibody

Specificity

Oryza sativa

Antibody description

Rabbit polyclonal antibody to LOC4331832

Preparation

This antigen of this antibody was recombinant protein within rice tify11a aa 50-179 / 179.

Formulation

Liquid, 1*TBS (pH7.4), 0.2% BSA, 50% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ELISA

Dilutions

WB: 1:500-1:2,000

ELISA: Use at an assay dependent concentration.

Validations




Fig1:; Western blot analysis of TIFY11A on rice tissue lysates. Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (1/1,000) was used in 5% NFDM/TBST at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1:200,000 dilution was used for 1 hour at room temperature.