

# SULT1A1 antibody



Catalog Number: 107612

## Product name

SULT1A1 antibody

## Specificity

Human; other species not tested.

## Antibody description

SULT1A1 Mouse Monoclonal antibody. Positive WB detected in Human placenta, HEK-293 cells, human liver tissue. Observed molecular weight by Western-blot: 34kd

## Preparation

This antibody was obtained by immunization of SULT1A1 recombinant protein (Accession Number: XM\_017023611). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:107612(SULT1A1 antibody) at dilution of 1:200