

Anti-FLI1 Rabbit antibody



Catalog Number: 168612

Product name

Anti-FLI1 Rabbit antibody

Specificity

Human, Mouse

Antibody description

FLI1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human FLI1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000

IHC: 1/50

Validations



Western blot detection of FLI1 in Jurkat, C6 cell lysates using FLI1 Rabbit mAb(1:500 diluted). Predicted band size: 51 kDa. Observed band size: 51 kDa.