

# Mouse Pcsk9/Narc1 (His tag) recombinant protein



Catalog Number: 518612

## General Information

### Gene Name Synonym

Proprotein Convertase Subtilisin/Kexin Type 9, Neural Apoptosis-Regulated Convertase 1, NARC-1, Proprotein Convertase 9, PC9, Subtilisin/Kexin-Like Protease PC9, Pcsk9, Narc1

### Protein Construction

A DNA sequence encoding the Mouse Pcsk9/Narc1 Gln35-Gln692 is expressed with a 6His tag at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 50mM HEPES, 150mM NaCl, pH 7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Mouse Pcsk9/Narc1 (His tag) recombinant protein