

Mouse Pcsk9/Narc1 (His tag) recombinant protein



Catalog Number: 518612

General Information

Gene Name Synonym

Proprotein Convertase Subtilisin/Kexin Type 9, Neural Apoptosis-Regulated Convertase 1, NARC-1, Proprotein Convertase 9, PC9, Subtilisin/Kexin-Like Protease PC9, Pcsk9, Narc1

Protein Construction

A DNA sequence encoding the Mouse Pcsk9/Narc1 Gln35-Gln692 is expressed with a 6His tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 50mM HEPES, 150mM NaCl, pH 7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Pcsk9/Narc1 (His tag) recombinant protein