

TRPM4 Polyclonal Antibody



Catalog Number: 160312

Product name

TRPM4 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to TRPM4

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human TRPM4 (NP_060106.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:200 - 1:1000

IHC 1:20 - 1:100

Validations



Western blot - TRPM4 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TRPM4 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 45s.