

# TRIM33 antibody



Catalog Number: 116312

---

## Product name

TRIM33 antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

Human; other species not tested.

## Applications

ELISA, WB

## Antibody description

TRIM33 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, MCF7 cells, PC-3 cells. Observed molecular weight by Western-blot: 140-150 kDa

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_015906). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:116312(TRIM33 antibody) at dilution of 1:500

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality