

# ATP6V1E2 antibody



Catalog Number: 108312

## Product name

ATP6V1E2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ATP6V1E2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, MCF7 cells. Observed molecular weight by Western-blot: 31-33kd

## Preparation

This antibody was obtained by immunization of ATP6V1E2 recombinant protein (Accession Number: XM\_011533153). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108312(ATP6V1E2 antibody) at dilution of 1:400