

# CIB4 antibody

Catalog Number: 109312

## Product name

CIB4 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CIB4 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 22 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001029881). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:109312(CIB4 antibody) at dilution of 1:100



Immunohistochemical of paraffin-embedded human kidney using Catalog No:109312(CIB4 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:109312(CIB4 antibody) at dilution of 1:50 (under 40x lens)