

Human PRND/DPL/UNQ1830/PRO3443 (His tag) recombinant protein



Catalog Number: 519712

General Information

Gene Name Synonym

Prion-like protein doppel, PrPLP, Prion protein 2, PRND, DPL

Protein Construction

A DNA sequence encoding the Human PRND/DPL/UNQ1830/PRO3443 Arg27-Gly152 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human PRND/DPL/UNQ1830/PRO3443 (His tag)
recombinant protein