

CD30 antibody

Catalog Number: 109022

Product name

CD30 antibody

Specificity

Human; other species not tested.

Antibody description

CD30 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue, human tonsillitis tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001243). Purification method: Antigen affinity purified.

Formulation

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human T-cell anaplastic large cell lymphoma (ALCL) tissue slide using Catalog No:109022(TNFRSF8 Antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human T-cell anaplastic large cell lymphoma (ALCL) tissue slide using Catalog No:109022(TNFRSF8 Antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieved with Tris-EDTA buffer(pH9).