

# EYA3 antibody

Catalog Number: 110422

## Product name

EYA3 antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

## Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

## Antibody description

EYA3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, L02 cells. Positive IP detected in L02 cells. Positive IF detected in HeLa cells. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 63-66 kDa

## Validations



## Preparation

This antibody was obtained by immunization of EYA3 recombinant protein (Accession Number: NM\_172098). Purification method: Antigen affinity purified.

Immunofluorescent analysis of HeLa cells, using EYA3 antibody Catalog No:110422 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, IF, IHC, WB, IP


## Dilutions

Recommended Dilution:

WB: 1:200-1:2000



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110422(EYA3 Antibody) at dilution of 1:600




# EYA3 antibody

Catalog Number: 110422


---

IP Result of anti-EYA3 (IP:Catalog No:110422, 4ug; Detection:Catalog No:110422 1:1000) with L02 cells lysate 1800ug.



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog

No:110422(EYA3 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:110422(EYA3 Antibody) at dilution of 1:50 (under 40x lens)