

Mouse LIFR/CD118 (His Tag) recombinant protein



Catalog Number: 501422

General Information

Gene Name Synonym

D-factor/LIF receptor

Protein Construction

A DNA sequence encoding the extracellular domain of mouse LIFR (NP_038612.1) (Met 1-Ser 828) was expressed, with a polyhistidine tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Activity

1. Measured by its ability to bind mouse LIF-Fc (Cat:50755-M02H) in a functional ELISA.
2. Measured by its ability to inhibit the recombinant human LIF mediated inhibition in the M1 mouse myeloid leukemia cells. The ED₅₀ for this effect is typically 5-20 ng/mL in the presence of 2 ng/mL recombinant human LIF.

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Leu 44

Molecular Mass

The recombinant mouse LIFR consists of 796 amino acids and has a predicted molecular mass of 90 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rm LIFR is approximately 110-120 kDa due to glycosylation.

Formulation

- Lyophilized from sterile PBS, pH 7.4
1. 5 % trehalose and mannitol are added as protectants before lyophilization.
 2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse LIFR / CD118 Protein (His Tag) SDS-PAGE