

HLA-DQA1 antibody



Catalog Number: 111422

Product name

HLA-DQA1 antibody

Specificity

Human; other species not tested.

Antibody description

HLA-DQA1 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of HLA-DQA1 recombinant protein (Accession Number: BC008585). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:111422(HLA-DQA1 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111422(HLA-DQA1 Antibody) at dilution of 1:200 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111422(HLA-DQA1 Antibody) at dilution of 1:200 (under 40x lens). heat mediated antigen

HLA-DQA1 antibody



Catalog Number: 111422

retrieved with Tris-EDTA buffer(pH9).