

GYPB Polyclonal Antibody



Catalog Number: 162422

Product name

GYPB Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to GYPB

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-91 of human GYPB (NP_002091.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - GYPB Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using GYPB antibody at dilution of 1:100. Blue: DAPI for nuclear staining.