

Anti-MACC1 antibody

Catalog Number: 176422

Product name

Anti-MACC1 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to MACC1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human macc1. aa range:411-460

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations



Fig1.; Western Blot analysis of L929 cells using MACC1 Polyclonal Antibody.. Secondary antibody [catalog#: HA1001) was diluted at 1:20000