

# Anti-MACC1 antibody



Catalog Number: 176422

## Product name

Polyclonal

Anti-MACC1 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse

## Applications

## Antibody description

WB

Rabbit polyclonal antibody to MACC1

## Dilutions

## Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human macc1. aa range:411-460

## Validations

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality



Fig1.; Western Blot analysis of L929 cells using MACC1 Polyclonal Antibody.. Secondary antibody [catalog#: HA1001) was diluted at 1:20000