

TRX2 antibody

Catalog Number: 116422

Product name

TRX2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TRX2 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in Raji cells, human liver tissue. Observed molecular weight by Western blot: 12 kDa

Preparation

This antibody was obtained by immunization of TRX2 recombinant protein (Accession Number: NM_012473). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:116422(TRX2, TXN2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:116422(TRX2, TXN2 antibody) at dilution of 1:100 (under 10x lens)