

Anti-Cdc25C Mouse antibody



Catalog Number: 168422

Product name

Anti-Cdc25C Mouse antibody

Specificity

Human

Antibody description

Cdc25C Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human CDC25C expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IHC-P


Dilutions

WB: 1/500 - 1/2000

IHC: 1/200 - 1/1000

ELISA: 1/10000

Validations



Western blot detection of Cdc25C in 293T, HeLa and K562 cell lysates using Cdc25C mouse mAb (1:500 diluted). Predicted band size: 60/75KDa. Observed band size: 75KDa.