

BCL2A1 antibody



Catalog Number: 108422

Product name

BCL2A1 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

BCL2A1 Rabbit Polyclonal antibody. Positive WB detected in human lung tissue. Observed molecular weight by Western-blot: 20 kDa

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of BCL2A1 recombinant protein (Accession Number: NM_004049). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:117095(BCL2A1 antibody) at dilution of 1:500

Clonality