

# CNDP2 antibody



Catalog Number: 109422

## Product name

CNDP2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CNDP2 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, human colon tissue, mouse testis tissue. Positive IF detected in PC-3 cells. Observed molecular weight by Western-blot: 54 kDa

## Preparation

This antibody was obtained by immunization of CNDP2 recombinant protein (Accession Number: XM\_005266728). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109422(CNDP2 antibody) at dilution of 1:200



Immunofluorescent analysis of PC-3 cells using Catalog No:109422(CNDP2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)