

# GADD45GIP1 antibody



Catalog Number: 110822

## Product name

GADD45GIP1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GADD45GIP1 Rabbit Polyclonal antibody. Positive IP detected in mouse heart tissue. Positive WB detected in mouse heart tissue, HEK-293 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 25-28kd

## Preparation

This antibody was obtained by immunization of GADD45GIP1 recombinant protein (Accession Number: NM\_052850). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110822(GADD45GIP1 antibody) at dilution of



1:300

IP Result of anti-GADD45GIP1 (IP:Catalog No:110822, 4ug; Detection:Catalog No:110822 1:400) with mouse heart tissue lysate 4000ug.