

# IQCD antibody

Catalog Number: 111822

---

**Product name**

IQCD antibody

**Specificity**

Human, Mouse; other species not tested.

**Antibody description**

IQCD Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 52 kDa

**Preparation**

This antibody was obtained by immunization of IQCD recombinant protein (Accession Number: XM\_047428263). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

ELISA, WB

**Dilutions**

Recommended Dilution:

WB: 1:500-1:5000

**Validations**

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:111822(IQCD Antibody) at dilution of 1:1000