

Anti-LOC102724594 antibody



Catalog Number: 175822

Product name

Anti-LOC102724594 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to LOC102724594

Preparation

This antigen of this antibody was synthesized peptide derived from human protein . at aa range: 60-140

Formulation

Liquid, Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.;; Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4°over night.