

Anti-GALT Rabbit antibody



Catalog Number: 167822

Product name

Anti-GALT Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

GALT Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human GALT

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1/1000

IP: 1/20

Validations



Western blot detection of GALT in HeLa, C6, Rat Brain, Ramos cell lysates using GALT Rabbit mAb(1:1000 diluted). Predicted band size:43kDa. Observed band size:43kDa.