

DUSP9 antibody



Catalog Number: 110122

Product name

DUSP9 antibody

Specificity

human, mouse; other species not tested.

Antibody description

DUSP9 Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive WB detected in mouse liver tissue, HepG2 cells. Observed molecular weight by Western-blot: 42-44 kDa

Preparation

This antibody was obtained by immunization of DUSP9 recombinant protein (Accession Number: XM_011531124). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110122(DUSP9 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using DUSP9 antibody Catalog No:110122 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).