

PIK3CD Polyclonal Antibody



Catalog Number: 163122

Product name

PIK3CD Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to PIK3CD

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 150-250 of human PIK3CD (NP_005017.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - PIK3CD Polyclonal Antibody

Western blot analysis of extracts of mouse thymus, using PIK3CD antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3min.