

Rhesus IL17BR / IL17RB / IL-17 Receptor B (Fc Tag) recombinant protein



Catalog Number: 504122

General Information

Protein Construction

A DNA sequence encoding the rhesus IL17RB (EHH16291.1) (Met1-Gly288) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Rhesus

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Arg 18

Molecular Mass

The recombinant rhesus IL17RB consists of 512 amino acids and predicts a molecular mass of 57

kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rhesus IL17BR / IL17RB / IL-17 Receptor B Protein (Fc Tag)