

# CCL8 Polyclonal Antibody



Catalog Number: 165122

## Product name

CCL8 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to CCL8

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 23-99 of human CCL8 (NP\_005614.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

## Validations



Immunohistochemistry - CCL8 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human liver using CCL8 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - CCL8 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate cancer using CCL8 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - CCL8 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using CCL8 antibody at dilution of

# CCL8 Polyclonal Antibody



Catalog Number: 165122

---

1:100 (40x lens).