

AP3B2 antibody

Catalog Number: 108122

Product name

AP3B2 antibody

Immunogen

[Human AP3B2 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AP3B2 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue, mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 120-140 kDa

Preparation

This antibody was obtained by immunization of AP3B2 recombinant protein (Accession Number: NM_004644). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108122(AP3B2 Antibody) at dilution of 1:600



IP Result of anti-AP3B2 (IP:Catalog No:108122, 4ug; Detection:Catalog No:108122 1:500) with mouse brain tissue lysate 3600ug.