

Human ZNHIT1/CGBP1/ZNFN4A1 (His tag) recombinant protein



Catalog Number: 518122

General Information

Gene Name Synonym

inc finger HIT domain-containing protein 1, cyclin-G1-binding protein 1, zinc finger protein subfamily 4A member 1, p18 hamlet, CGBP1, ZNFN4A1, ZNHIT1

Protein Construction

A DNA sequence encoding the Human ZNHIT1/CGBP1/ZNFN4A1 Met1-Val154 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB \square 150mM NaCl \square 50% Glycerol,pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human ZNHIT1/CGBP1/ZNFN4A1 (His tag)
recombinant protein