

# Human ZNHIT1/CGBP1/ZNFN4A1 (His tag) recombinant protein



Catalog Number: 518122

## General Information

### Gene Name Synonym

inc finger HIT domain-containing protein 1, cyclin-G1-binding protein 1, zinc finger protein subfamily 4A member 1, p18 hamlet, CGBP1, ZNFN4A1, ZNHIT1

### Protein Construction

A DNA sequence encoding the Human ZNHIT1/CGBP1/ZNFN4A1 Met1-Val154 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB□ 150mM NaCl□50% Glycerol,pH7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human ZNHIT1/CGBP1/ZNFN4A1 (His tag)  
recombinant protein