

Anti-Cyclin E1 Mouse antibody



Catalog Number: 169122

Product name

Anti-Cyclin E1 Mouse antibody

Specificity

Human

Antibody description

Cyclin E1 Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human CCNE1 (AA: 307-410) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

IHC-P, FC


Validations



Immunohistochemical analysis of paraffin-embedded rectum cancer tissues using CCNE1 mouse mAb with DAB staining.



Flow cytometric analysis of K562 cells using CCNE1 mouse mAb (green) and negative control



Western blot analysis using CCNE1 mAb against

Anti-Cyclin E1 Mouse antibody



Catalog Number: 169122

Western blot analysis using CCNE1 mAb against HEK293 (1) and CCNE1 (AA)

