

# Anti-Cyclin E1 Mouse antibody



Catalog Number: 169122

## Product name

Anti-Cyclin E1 Mouse antibody

## Specificity

Human

## Antibody description

Cyclin E1 Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human CCNE1 (AA: 307-410) expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Mouse IgG1


## Applications

IHC-P, FC


## Validations



Immunohistochemical analysis of paraffin-embedded rectum cancer tissues using CCNE1 mouse mAb with DAB staining.



Flow cytometric analysis of K562 cells using CCNE1 mouse mAb (green) and negative control



Western blot analysis using CCNE1 mAb against

# Anti-Cyclin E1 Mouse antibody



Catalog Number: 169122

---

Western blot analysis using CCNE1 mAb against  
HEK293 (1) and CCNE1 (AA

