

METTL7B antibody



Catalog Number: 112522

Product name

METTL7B antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

METTL7B Rabbit Polyclonal antibody. Positive IHC detected in human liver. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 28kd

Preparation

This antibody was obtained by immunization of METTL7B recombinant protein (Accession Number: BC020509). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:112522(METTL7B antibody) at dilution of



Immunohistochemical of paraffin-embedded human liver using Catalog No:112522(METTL7B antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:112522(METTL7B antibody) at dilution of 1:100 (under 40x lens)