

# OVOL1 antibody



Catalog Number: 113522

## Product name

OVOL1 antibody

## Immunogen

[Human OVOL1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

OVOL1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse pancreas tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 30kd

## Preparation

This antibody was obtained by immunization of OVOL1 recombinant protein (Accession Number: NM\_004561). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113522(OVOL1 antibody) at dilution of 1:400



IP Result of anti-OVOL1 (IP:Catalog No:113522, 4ug; Detection:Catalog No:113522 1:500) with mouse brain tissue lysate 4000ug.