

Mouse BACE1 (His Tag) recombinant protein



Catalog Number: 504522

General Information

Protein Construction

A DNA sequence encoding the extracellular domain (Met 1-Thr 457) of mouse BACE1 (NM_011792.5) precursor was expressed with a C-terminal polyhistidine tag.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Activity

Measured by its ability to cleave a fluorescent peptide substrate Mca-Ser-Glu-Val-Asn-Leu-Asp-Ala-Glu-Phe-Arg-Lys(Dpn)-Arg-Arg-NH₂ (Catalog# ES004, R&D Systems). Cleavage of ES004 can be measured using excitation and emission wavelengths of 320 and 405 nm, respectively. The specific activity is >2 pmoles/min/ μ g.

Purity

> 97 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Thr 22

Molecular Mass

The secreted recombinant mouse BACE1 consists of 447 amino acids and has a calculated molecular mass of 49.8 kDa. As a result of glycosylation, the recombinant protein migrates as an approximately 60-65 kDa protein in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse BACE-1 Protein (His Tag) SDS-PAGE