

Anti-CHURC1-FNTB antibody



Catalog Number: 176522

Product name

Anti-CHURC1-FNTB antibody

Specificity

Human, Rat

Antibody description

Rabbit monoclonal antibody to CHURC1-FNTB

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*TBS (pH7.4), 0.05% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF

Dilutions

WB: 1:1,000-1:2,000

ICC/IF: 1:50-1:200

Validations





Fig1: Western blot analysis of FNTB on K562 cells lysates using anti-FNTB antibody at 1/1,000



dilution.

Fig2: ICC staining FNTB in JAR cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.