

Anti-HDAC3 (7G6) Mouse antibody



Catalog Number: 168522

Product name

Anti-HDAC3 (7G6) Mouse antibody

Specificity

Human

Antibody description

HDAC3 (7G6) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human HDAC3 (aa224-428) expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

WB, ICC/IF, IP


Dilutions

WB: 1/500 - 1/2000


ICC: 1/200 - 1/1000

ELISA: 1/10000


Validations



Immunoprecipitation analysis of HeLa cell lysates using HDAC3 mouse mAb.



Western blot detection of HDAC3 in COS7, C6, HeLa and 3T3 cell lysates using HDAC3 mouse mAb (dilution 1:500). Predicted band size: 48.8 kDa. Observed band size: 48.8 kDa.



Immunocytochemistry of HeLa cells using HDAC3 mouse mAb (dilution 1:100).