

KCNA10 antibody

Catalog Number: 111922

Product name

KCNA10 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KCNA10 Rabbit Polyclonal antibody. Positive IHC detected in human placenta tissue.

Preparation

This antibody was obtained by immunization of KCNA10 recombinant protein (Accession Number: NM_005549). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog No:111922(KCNA10 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog No:111922(KCNA10 Antibody) at dilution of 1:50 (under 40x lens)