

GJC2 Polyclonal Antibody



Catalog Number: 163922

Product name

GJC2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to GJC2

Preparation

Antigen: A synthetic peptide of human GJC2

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - GJC2 Polyclonal Antibody

Western blot analysis of extracts of mouse brain, using GJC2 antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time:



30s.

Immunohistochemistry - GJC2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human kidney cancer using GJC2 antibody at dilution of 1:200 (40x lens).