

Anti-Arylsulfatase A / ARSA antibody



Catalog Number: 103922

Product name

Anti-Arylsulfatase A / ARSA antibody

Immunogen

[Human Arylsulfatase A / ARSA \(His Tag\) recombinant protein](#)

Specificity

Human Arylsulfatase A / ARSA

Has cross-reactivity in ELISA with Mouse ARSA

Antibody description

Rabbit monoclonal to Arylsulfatase A / ARSA

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human Arylsulfatase A / ARSA (rh Arylsulfatase A / ARSA; NP_000478.2; Met 1-Ala 507).

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 2-10 µg/mL

Validations



Items	Lanes	A
Sample (whole cell lysate)		HepG2
Sample Volume (µg/lane)		60
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		2-10 µg/ml
Secondary Antibody		Dylight 800 labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.

Arylsulfatase A / ARSA Antibody, Rabbit MAb, Western blot