

# DHRS11 antibody

Catalog Number: 109922

## Product name

DHRS11 antibody

## Immunogen

[Human DHRS11 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DHRS11 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, COLO 320 cells. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 25-28kd

## Preparation

This antibody was obtained by immunization of DHRS11 recombinant protein (Accession Number: NM\_024308). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109922(DHRS11 antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:109922(DHRS11 Antibody) at dilution of 1:50 (under 40x lens)